

Product datasheet for **TL519685**

Rp1l1 Mouse shRNA Plasmid (Locus ID 271209)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Rp1l1 Mouse shRNA Plasmid (Locus ID 271209)
Locus ID:	271209
Synonyms:	BC031365; BE983540; Dcdc4; Rp1; Rp1hl1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Rp1l1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 271209). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC129887 , NM_146246 , NM_146246.2 , NM_146246.3 , BC031365 , BC049924
UniProt ID:	Q8CGM2
Summary:	This gene encodes a member of the doublecortin family. This protein is a retinal-specific protein. It contains two N-terminal doublecortin domains, which can assemble and stabilize axonemal microtubules, but lacks the C-terminal repetitive regions including the 16aa repeat found in human RP1L1. This protein and the RP1 protein, another retinal-specific protein, play essential and synergistic roles in affecting photosensitivity and outer segment morphogenesis of rod photoreceptors. [provided by RefSeq, Sep 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).