

Product datasheet for **TL701809V**

Akr1c2 Rat shRNA Lentiviral Particle (Locus ID 291283)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Akr1c2 Rat shRNA Lentiviral Particle (Locus ID 291283)
Locus ID:	291283
Synonyms:	Akr1c21
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Akr1c21 - Rat shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_001013057 , NM_001013057.1 , BC078957
UniProt ID:	Q6AYQ2
Summary:	NADP-dependent 17-alpha-hydroxysteroid dehydrogenase that converts 5-alpha-androstane-3,17-dione into androsterone. Has lower 3-alpha-hydroxysteroid dehydrogenase activity. Has broad substrate specificity and acts on various 17-alpha-hydroxysteroids, 17-ketosteroids, 3-alpha hydroxysteroids and 3-ketosteroids. Reduction of keto groups is strictly stereoselective. Reduction of 17-ketosteroids yields only 17-alpha-hydroxysteroids. Likewise, reduction of 3-ketosteroids yields only 3-alpha-hydroxysteroids (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).