

## Product datasheet for **TL701987V**

### Lsg1 Rat shRNA Lentiviral Particle (Locus ID 288029)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Lsg1 Rat shRNA Lentiviral Particle (Locus ID 288029)
Locus ID:	288029
Synonyms:	RGD1309089
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Lsg1 - Rat shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_001013421</a> , <a href="#">NM_001013421.1</a> , <a href="#">BC079013</a> , <a href="#">BC091338</a>
UniProt ID:	<a href="#">Q5BJT6</a>
Summary:	GTPase required for the XPO1/CRM1-mediated nuclear export of the 60S ribosomal subunit. Probably acts by mediating the release of NMD3 from the 60S ribosomal subunit after export into the cytoplasm (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .
Performance Guaranteed:	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at <a href="mailto:techsupport@origene.com">techsupport@origene.com</a>. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>



[View online »](#)