

## Product datasheet for **TL703963**

### Shtn1 Rat shRNA Plasmid (Locus ID 292139)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Shtn1 Rat shRNA Plasmid (Locus ID 292139)
Locus ID:	292139
Synonyms:	RGD1311558; Shootin1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Shtn1 - Rat, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 292139). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001079705</a> , <a href="#">NM_001303537</a> , <a href="#">NM_001079705.1</a> , <a href="#">NM_001079705.2</a> , <a href="#">NM_001079705.4</a> , <a href="#">BC085922</a>
UniProt ID:	<a href="#">A0MZ67</a>
Summary:	Involved in the generation of internal asymmetric signals required for neuronal polarization and neurite outgrowth (PubMed:17030985, PubMed:17439943, PubMed:18519736, PubMed:20664640). Mediates netrin-1-induced F-actin-substrate coupling or 'clutch engagement' within the axon growth cone through activation of CDC42, RAC1 and PAK1-dependent signaling pathway, thereby converting the F-actin retrograde flow into traction forces, concomitantly with filopodium extension and axon outgrowth (PubMed:18519736, PubMed:23453953). Plays a role in cytoskeletal organization by regulating the subcellular localization of phosphoinositide 3-kinase (PI3K) activity at the axonal growth cone (PubMed:17030985). Plays also a role in regenerative neurite outgrowth (PubMed:20664640). In the developing cortex, cooperates with KIF20B to promote both the transition from the multipolar to the bipolar stage and the radial migration of cortical neurons from the ventricular zone toward the superficial layer of the neocortex. Involved in the accumulation of phosphatidylinositol 3,4,5-trisphosphate (PIP3) in the growth cone of primary hippocampal neurons (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact [techsupport@origene.com](mailto:techsupport@origene.com). If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).