

## Product datasheet for **TL703977**

### Sephs2 Rat shRNA Plasmid (Locus ID 308993)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Sephs2 Rat shRNA Plasmid (Locus ID 308993)
Locus ID:	308993
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Sephs2 - Rat, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 308993). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001079889</a> , <a href="#">NM_001079889.1</a> , <a href="#">NM_001079889.2</a> , <a href="#">BC081915</a> , <a href="#">BC097429</a> , <a href="#">BC128783</a>
Summary:	This gene encodes an enzyme that catalyzes the production of monoselenophosphate (MSP) from selenide and ATP. MSP is the selenium donor required for synthesis of selenocysteine (Sec), which is co-translationally incorporated into selenoproteins at in-frame UGA codons that normally signal translation termination. The 3' UTRs of selenoprotein mRNAs contain a conserved stem-loop structure, the Sec insertion sequence (SECIS) element, which is necessary for the recognition of UGA as a Sec codon rather than as a stop signal. This protein is itself a selenoprotein containing a Sec residue at its active site, suggesting the existence of an autoregulatory mechanism. It is preferentially expressed in tissues implicated in the synthesis of selenoproteins and in sites of blood cell development. [provided by RefSeq, May 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).