

## Product datasheet for **TL709756V**

### Cckbr Rat shRNA Lentiviral Particle (Locus ID 25706)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Cckbr Rat shRNA Lentiviral Particle (Locus ID 25706)
Locus ID:	25706
Synonyms:	Cck2r; Cholrec
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Cckbr - Rat shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_013165</a> , <a href="#">NM_013165.1</a> , <a href="#">NM_013165.2</a>
UniProt ID:	<a href="#">P30553</a>
Summary:	Cholecystokinin was one of the first gastrointestinal peptides discovered in the mammalian brain. Cholecystokinin receptors are members of the G protein-coupled receptor superfamily, stimulating phosphatidylinositol turnover and intracellular calcium mobilization. This gene encodes a cholecystokinin receptor of the B subtype. These receptors occur throughout the central nervous system where they modulate anxiety, analgesia, arousal, and neuroleptic activity. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).