

Product datasheet for **TL710112V**

Mapk9 Rat shRNA Lentiviral Particle (Locus ID 50658)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Mapk9 Rat shRNA Lentiviral Particle (Locus ID 50658)
Locus ID:	50658
Synonyms:	SAPK
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Mapk9 - Rat shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_001270544 , NM_001270545 , NM_017322 , NM_017322.1 , NM_017322.2 , NM_001270545.1 , NM_001270544.1 , BC061870
UniProt ID:	P49186
Summary:	The protein encoded by this gene is a member of the MAP kinase family. MAP kinases are multifunctional proteins. They are involved in a wide variety of cellular processes such as growth, proliferation, differentiation, transcription regulation, and development. They function in stress responses, apoptosis, inflammation and transformation. For example, this kinase plays a dominant role in mediating proliferation of lung cancer and prostate cancer cells, and in regulation of osmotic stress-induced tight junction disruption. Several alternatively spliced transcript variants encoding distinct isoforms have been reported. [provided by RefSeq, Jul 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).