

Product datasheet for TL710189

Snca Rat shRNA Plasmid (Locus ID 29219)

Product data:

Product Type: shRNA Plasmids

Product Name: Snca Rat shRNA Plasmid (Locus ID 29219)

Locus ID: 29219

Synonyms: MGC105443

Vector: pGFP-C-shLenti (TR30023)

E. coli Selection: Chloramphenicol (34 ug/ml)

Mammalian Cell

Selection:

Puromycin

Format: Lentiviral plasmids

Components: Snca - Rat, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 29219). 5µg

purified plasmid DNA per construct

29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.

RefSeq: NM 019169, NM 019169.1, NM 019169.2, BC087682

UniProt ID: P37377

Summary: may play a role in regulation of synaptic vesicle biogenesis and organization; may be involved

in neurotransmission [RGD, Feb 2006]

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To

be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u>. If you need a special design or shRNA sequence, please utilize our <u>custom shRNA service</u>.

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com




Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:




GFP signal was observed under microscope at 48 hours after transduction of TL710189A virus into HEK293 cells. TL710189A virus was prepared using lenti-shRNA TL710189A and [TR30037] packaging kit.






GFP signal was observed under microscope at 48 hours after transduction of TL710189B virus into HEK293 cells. TL710189B virus was prepared using lenti-shRNA TL710189B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL710189C] virus into HEK293 cells. [TL710189C] virus was prepared using lenti-shRNA [TL710189C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL710189D] virus into HEK293 cells. [TL710189D] virus was prepared using lenti-shRNA [TL710189D] and [TR30037] packaging kit.