

Product datasheet for **TL710880**

Chd8 Rat shRNA Plasmid (Locus ID 65027)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Chd8 Rat shRNA Plasmid (Locus ID 65027)
Locus ID:	65027
Synonyms:	CHD-8
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Chd8 - Rat, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 65027). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001347661 , NM_022933 , NM_022933.1 , NM_022933.2
Summary:	This gene encodes a member of the chromodomain-helicase-DNA binding protein family, which is characterized by a SNF2-like domain and two chromatin organization modifier domains. The encoded protein also contains brahma and kismet domains, which is common to the subfamily of chromodomain-helicase-DNA binding proteins to which this protein belongs. In mammals, this gene has been shown to function in several processes including transcriptional regulation, epigenetic remodeling, promotion of cell proliferation, and regulation of RNA synthesis. Knockout of this gene in mice causes early embryonic lethality due to widespread apoptosis. Heterozygous loss of function mutations in mice result in autism spectrum disorder-like behaviors that include increased anxiety, repetitive behavior, and altered social behavior. [provided by RefSeq, Dec 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).