

Product datasheet for **TL711914**

Hrh3 Rat shRNA Plasmid (Locus ID 85268)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Hrh3 Rat shRNA Plasmid (Locus ID 85268)
Locus ID:	85268
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Hrh3 - Rat, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 85268). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001270564 , NM_001270565 , NM_001270566 , NM_001270567 , NM_001270568 , NM_001270569 , NM_001270570 , NM_053506 , NR_073042 , NM_053506.1 , NM_053506.2 , NM_001270570.1 , NM_001270565.1 , NM_001270569.1 , NM_001270564.1 , NM_001270568.1 , NM_001270567.1 , NM_001270566.1 , BC087707
UniProt ID:	Q9QYN8
Summary:	This gene encodes a histamine H3 receptor that belongs to the superfamily of G-protein coupled receptors. This protein functions as a presynaptic autoreceptor on histamine neurons in the brain, and a presynaptic heteroreceptor in nonhistamine-containing neurons in both the central and peripheral nervous systems. It is deemed a great target for the development of therapeutics for numerous disorders, including obesity, epilepsy, and such cognitive diseases as attention deficit hyperactivity disorder and Alzheimer's disease. Several alternatively spliced transcript variants encoding different isoforms, with different brain expression patterns and signaling properties, have been described for this gene. [provided by RefSeq, Jul 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).