

Product datasheet for **TL712174V**

Nrxn3 Rat shRNA Lentiviral Particle (Locus ID 116508)


Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Nrxn3 Rat shRNA Lentiviral Particle (Locus ID 116508)
Locus ID:	116508
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Nrxn3 - Rat shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_053817 , NM_053817.1 , NM_053817.2
UniProt ID:	Q07310
Summary:	brain specific cell surface protein; may act as a receptor and play a role in cellular signalling [RGD, Feb 2006]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .
Performance Guaranteed:	OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).




[View online »](#)

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL712174A virus into HEK293 cells. TL712174A virus was prepared using lenti-shRNA TL712174A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL712174B virus into HEK293 cells. TL712174B virus was prepared using lenti-shRNA TL712174B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL712174C] virus into HEK293 cells. [TL712174C] virus was prepared using lenti-shRNA [TL712174C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL712174D] virus into HEK293 cells. [TL712174D] virus was prepared using lenti-shRNA [TL712174D] and [TR30037] packaging kit.