

Product datasheet for TL900004V

KD-Validated (via WB) CDKN1B Lentiviral shRNA Knockdown Kit

Product data:

Product Type: Lenti shRNA Knockdown Kits

Tested Cell Line: HSHC

Components: 3 1.5mL vials

1. WB-validated CDKN1B shRNA lentiviral particles (0.5mL, $> 1 \times 10^7$ TU/mL)
2. Scrambled negative control shRNA lentiviral particles (0.5mL, $> 1 \times 10^7$ TU/mL)
3. Verification Tool: KD-Validated Anti-Cyclin Dependent Kinase Inhibitor 1B Rabbit mAb TA424702 (5 μ L)

Storage: The lentiviral particles should be stored at -80°C . The antibody should be stored at -20°C .

Product images:

