

Product datasheet for TL900012V

KD-Validated (via WB) UCHL1 Lentiviral shRNA Knockdown Kit

Product data:

Product Type: Lenti shRNA Knockdown Kits


Tested Cell Line: HSHC

Components: 3 1.5mL vials

1. WB-validated UCHL1 shRNA lentiviral particles (0.5mL, $> 1 \times 10^7$ TU/mL)
2. Scrambled negative control shRNA lentiviral particles (0.5mL, $> 1 \times 10^7$ TU/mL)
3. Verification Tool: KD-Validated Anti-UCHL1 Rabbit mAb TA425322 (5 μ L)

Storage: The lentiviral particles should be stored at -80°C The antibody should be stored at -20°C

Product images:



Western blotting analysis. UCHL1 protein expression in wild-type (WT) and shRNA knockdown (KD) HSHC cells was detected using Western blotting. β -Tubulin served as a loading control. The blots were incubated with primary antibodies against UCHL1 and β -Tubulin, respectively, followed by incubating with HRP-conjugated goat anti-rabbit secondary antibody. Images were developed using FeQ™ ECL Substrate Kit.

