

Product datasheet for TL900017V

KD-Validated (via WB) PUS7 Lentiviral shRNA Knockdown Kit

Product data:


Product Type: Lenti shRNA Knockdown Kits

Tested Cell Line: HeLa

Components: 3 1.5mL vials
1. WB-validated PUS7 shRNA lentiviral particles (0.5mL, $> 1 \times 10^7$ TU/mL)
2. Scrambled negative control shRNA lentiviral particles (0.5mL, $> 1 \times 10^7$ TU/mL)
3. Verification Tool: KD-Validated Anti-PUS7 Mouse mAb (5 μ L)

Storage: The lentiviral particles should be stored at -80°C . The antibody should be stored at -20°C .

Product images:



Western blotting analysis. PUS7 protein expression in wild-type (WT) and shRNA knockdown (KD) HeLa cells was detected using Western blotting. Hsp90 α served as a loading control. The blots were incubated with primary antibodies against PUS7 and Hsp90 α , respectively, followed by incubating with HRP-conjugated goat anti-mouse secondary antibody. Images were developed using NaQ™ ECL Substrate Kit.

