

Product datasheet for TL900023V

KD-Validated (via WB) FLOT2 Lentiviral shRNA Knockdown Kit

Product data:

Product Type: Lenti shRNA Knockdown Kits

Tested Cell Line: HepG2

Components: 3 1.5mL vials
1. WB-validated FLOT2 shRNA lentiviral particles (0.5mL, $> 1 \times 10^7$ TU/mL)
2. Scrambled negative control shRNA lentiviral particles (0.5mL, $> 1 \times 10^7$ TU/mL)
3. Verification Tool: KD-Validated Anti-Flotillin 2 Rabbit mAb (5 μ L)

Storage: The lentiviral particles should be stored at -80°C . The antibody should be stored at -20°C .

Product images:

