

Product datasheet for **TL900039V**

KD-Validated (via WB) RAF1 Lentiviral shRNA Knockdown Kit

Product data:

Product Type: Lenti shRNA Knockdown Kits


Tested Cell Line: HeLa

Components: 3 1.5mL vials

1. WB-validated RAF1 shRNA lentiviral particles (0.5mL, > 1×10⁷ TU/mL)
2. Scrambled negative control shRNA lentiviral particles (0.5mL, > 1×10⁷ TU/mL)
3. Verification Tool: KD-Validated Anti-Raf1 Rabbit mAb TA425184 (5 µL)

Storage: The lentiviral particles should be stored at -80°C. The antibody should be stored at -20°C.

Product images:



Western blotting analysis. RAF1 protein expression in wild-type (WT) and shRNA knockdown (KD) HeLa cells was detected using Western blotting. Hsp90 α served as a loading control. The blots were incubated with primary antibodies against RAF1 and Hsp90 α, respectively, followed by incubating with HRP-conjugated goat anti-rabbit secondary antibody. Images were developed using FeQ™ ECL Substrate Kit (Cat#226).



[View online »](#)