

Product datasheet for TL900049V

KD-Validated (via WB) IFT88 Lentiviral shRNA Knockdown Kit

Product data:

Product Type: Lenti shRNA Knockdown Kits

Tested Cell Line: HeLa

Components: 3 1.5mL vials


1. WB-validated IFT88 shRNA lentiviral particles (0.5mL, > 1×10⁷ TU/mL)

2. Scrambled negative control shRNA lentiviral particles (0.5mL, > 1×10⁷ TU/mL)

3. Verification Tool: KD-Validated Anti-Intraflagellar Transport 88 Rabbit mAb TA424899 (5 µL)

Storage: The lentiviral particles should be stored at -80°C. The antibody should be stored at -20°C.

Product images:



Western blotting analysis. IFT88 protein expression in wild-type (WT) and shRNA knockdown (KD) HeLa cells was detected using Western blotting. β-Tubulin served as a loading control. The blots were incubated with primary antibodies against IFT88 and β-Tubulin, respectively, followed by incubating with HRP-conjugated goat anti-rabbit secondary antibody. Images were developed using FeQ™ ECL Substrate Kit (Cat#226).



View online »