

Product datasheet for TL900068V

KD-Validated (via WB) MEF2A Lentiviral shRNA Knockdown Kit

Product data:

Product Type: Lenti shRNA Knockdown Kits


Tested Cell Line: HepG2

Components: 3 1.5mL vials

1. WB-validated MEF2A shRNA lentiviral particles (0.5mL, $> 1 \times 10^7$ TU/mL)
2. Scrambled negative control shRNA lentiviral particles (0.5mL, $> 1 \times 10^7$ TU/mL)
3. Verification Tool: KD-Validated Anti-Myocyte Enhancer Factor 2A Rabbit mAb (5 μ L)

Storage: The lentiviral particles should be stored at -80°C The antibody should be stored at -20°C

Product images:



Western blotting analysis. MEF2A protein expression in wild-type (WT) and shRNA knockdown (KD) HepG2 cells was detected using Western blotting. Hsp90 α served as a loading control. The blots were incubated with primary antibodies against MEF2A and Hsp90 α , respectively, followed by incubating with HRP-conjugated goat anti-rabbit secondary antibody. Images were developed using FeQ™ ECL Substrate Kit.



View online »