

Product datasheet for **TL900072V**

KD-Validated (via WB) FZD9 Lentiviral shRNA Knockdown Kit

Product data:

Product Type: Lenti shRNA Knockdown Kits


Tested Cell Line: HeLa

Components: 3 1.5mL vials

1. WB-validated FZD9 shRNA lentiviral particles (0.5mL, $> 1 \times 10^7$ TU/mL)
2. Scrambled negative control shRNA lentiviral particles (0.5mL, $> 1 \times 10^7$ TU/mL)
3. Verification Tool: KD-Validated Anti-Frizzled Class Receptor 9 Rabbit mAb TA424788 (5 μ L)

Storage: The lentiviral particles should be stored at -80°C . The antibody should be stored at -20°C .

Product images:



Western blotting analysis. FZD9 protein expression in wild-type (WT) and shRNA knockdown (KD) HeLa cells was detected using Western blotting. Hsp90 α served as a loading control. The blots were incubated with primary antibodies against FZD9 and Hsp90 α , respectively, followed by incubating with HRP-conjugated goat anti-rabbit secondary antibody. Images were developed using FeQ™ ECL Substrate Kit (Cat#226).



View online »