

Product datasheet for TL900073V

KD-Validated (via WB) BIRC2 Lentiviral shRNA Knockdown Kit

Product data:

Product Type: Lenti shRNA Knockdown Kits

Tested Cell Line: HSHC

Components: 3 1.5mL vials


1. WB-validated BIRC2 shRNA lentiviral particles (0.5mL, > 1×10⁷ TU/mL)

2. Scrambled negative control shRNA lentiviral particles (0.5mL, > 1×10⁷ TU/mL)

3. Verification Tool: KD-Validated Baculoviral IAP Repeat Containing 2 Recombinant Rabbit mAb (5 µL)

Storage: The lentiviral particles should be stored at -80°C. The antibody should be stored at -20°C.

Product images:



View online »