

## Product datasheet for TP300035

### EXOSC3 (NM\_016042) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human exosome component 3 (EXOSC3), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200035 protein sequence Red=Cloning site Green=Tags(s)
	<p>MAEPASVAAESLAGSRARAARTVLGQVVLPGEEELLPEQEDAEGPGGAVERPLSLNARACSRVRVWCGPG LRRCGDRLLVTKCGRLRHKEPGSGGGVYVWDSQQKRYVPVKGDHVGIVTAKSGDIFKVDVGGSEPAS LSYLSFEGATKRNRPNVQVGDLIYGFVWANKDMEPEMVCIDSCGRANGMGVIGQDGLLFKVTGLIRKL LAPDCEIIQEVGKLYPLEIVFGMNGRIWVKAKTIQQLILANILEACEHMTSDQRKQIFSRLAES</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	29.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_057126</a>
Locus ID:	51010
UniProt ID:	<a href="#">Q9NOT5</a>



[View online »](#)

RefSeq Size: 1857

Cytogenetics: 9p13.2

RefSeq ORF: 825


Synonyms: bA3J10.7; CGI-102; hRrp-40; p10; PCH1B; RRP40; Rrp40p

**Summary:** This gene encodes a non-catalytic component of the human exosome, a complex with 3'-5' exoribonuclease activity that plays a role in numerous RNA processing and degradation activities. Related pseudogenes of this gene are found on chromosome 19 and 21. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jun 2012]

**Protein Families:** Stem cell - Pluripotency

**Protein Pathways:** RNA degradation

### Product images:



Coomassie blue staining of purified EXOSC3 protein (Cat# TP300035). The protein was produced from HEK293T cells transfected with EXOSC3 cDNA clone (Cat# [RC200035]) using MegaTran 2.0 (Cat# [TT210002]).