

Product datasheet for **TP300036L**

HSPC152 (TRMT112) (NM_016404) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hypothetical protein HSPC152 (HSPC152), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200036 protein sequence Red =Cloning site Green =Tags(s)
	MKLLTHNLLSSHVRGVGSRGFPLRLQATEVTRICPVEFNPNFVARMIPKVEWSAFLEAADNLRLLIQVPGKP VEGYEENEEFLRTMHHLLEVEVIEGTLQCPEGRMFPISRGIPNMLLSEETES
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	14 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_057488
Locus ID:	51504
UniProt ID:	Q9UI30
RefSeq Size:	1305
Cytogenetics:	11q13.1




[View online »](#)

RefSeq ORF: 375

Synonyms: HSPC152; HSPC170; hTrm112; TRM112; TRMT11-2

Summary: Acts as an activator of both rRNA/tRNA and protein methyltransferases (PubMed:25851604). Together with methyltransferase BUD23, methylates the N(7) position of a guanine in 18S rRNA (PubMed:25851604). The heterodimer with HEMK2/N6AMT1 catalyzes N5-methylation of ETF1 on 'Gln-185', using S-adenosyl L-methionine as methyl donor (PubMed:18539146). The heterodimer with ALKBH8 catalyzes the methylation of 5-carboxymethyl uridine to 5-methylcarboxymethyl uridine at the wobble position of the anticodon loop in target tRNA species (PubMed:20308323). Involved in the pre-rRNA processing steps leading to small-subunit rRNA production (PubMed:25851604).[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified TRMT112 protein (Cat# [TP300036]). The protein was produced from HEK293T cells transfected with TRMT112 cDNA clone (Cat# [RC200036]) using MegaTran 2.0 (Cat# [TT210002]).