

Name: Recombinant protein of human ribosomal protein L36 (RPL36), transcript variant 2

Catalog: TP300054

Product Data Sheet

Gene Name: Homo sapiens ribosomal protein L36 (RPL36), transcript variant 2
GenBank accession: NM_015414

Description: Recombinant protein of human ribosomal protein L36 (RPL36), transcript variant 2

Protein Gel Analysis Data:

Amount: 20 ug

Buffer and Storage: 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol.

Source: Recombinant protein was produced with TrueORF clone, [RC200054](#), encoding the full-length human RPL36 with C-terminal DDK tag, from human HEK293 cells.

Protein Accn: NP_056229


Gene Synonym: L36; L36

Predicted molecular weight: 12.1 kDa

Tags: C-terminal MYC/DDK

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining.

Concentration: >50 ug/mL as determined by microplate BCA method



Purified recombinant protein RPL36 was analyzed by SDS-PAGE gel and Coomassie Blue Staining.

This product is for research and development use only. Not for use in human.