

## Product datasheet for **TP300075L**

### CUTA (NM\_001014838) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human cutA divalent cation tolerance homolog (E. coli) (CUTA), transcript variant 4, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200075 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MPALLPVASRLLLLPRVLLTMASGSPPTQPSPASDSGSGYVPGSVSAAFVTCPNEKVAKEIARAVVEKRL AACVNLIPQITSIYEWKKGKIEEDSEVLMMIKTQSSLVPALTDVRSVHPYEVAEVIALPVEQGNFPYLQW VRQVTESVSDSITVLP
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	16.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001014838</a>
Locus ID:	51596
UniProt ID:	<a href="#">O60888</a>



[View online »](#)

RefSeq Size:	1041
Cytogenetics:	6p21.32
RefSeq ORF:	468
Synonyms:	ACHAP; C6orf82
Summary:	May form part of a complex of membrane proteins attached to acetylcholinesterase (AChE). [UniProtKB/Swiss-Prot Function]
Protein Families:	Transmembrane

### Product images:



Coomassie blue staining of purified CUTA protein (Cat# [TP300075]). The protein was produced from HEK293T cells transfected with CUTA cDNA clone (Cat# [RC200075]) using MegaTran 2.0 (Cat# [TT210002]).