

Product datasheet for **TP300076L**

DBNDD2 (NM_001048221) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human dysbindin (dystrobrevin binding protein 1) domain containing 2 (DBNDD2), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200076 protein sequence Red =Cloning site Green =Tags(s)
	<p>MDPNPRAALERQQLRLRERQKFFEDILQPETEFVFPPLSHLHLESQRPPIGSISSEMEVNVDTLEQVELIDL GDPDAADVFLPCEDPPPTQSSGVDNHLEELSLPVPTSDRTTSRTSSSSSSSDSSTNLHSPNPSSDDGADTP LAQSDEEEERGDGGAEPGACS</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	17.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001041686
Locus ID:	55861
UniProt ID:	Q9BOY9 , Q9BOY9-2



[View online »](#)

RefSeq Size:	1165
Cytogenetics:	20q13.12
RefSeq ORF:	483
Synonyms:	C20orf35; CK1BP; HSMNP1
Summary:	May modulate the activity of casein kinase-1. Inhibits CSNK1D autophosphorylation (in vitro). [UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified DBNDD2 protein (Cat# [TP300076]). The protein was produced from HEK293T cells transfected with DBNDD2 cDNA clone (Cat# [RC200076]) using MegaTran 2.0 (Cat# [TT210002]).