

Product datasheet for TP300092

SKIP (INPP5K) (NM_016532) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human inositol polyphosphate-5-phosphatase K (INPP5K), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200092 protein sequence Red =Cloning site Green =Tags(s)
	<p>MSSRKLSGPKGRRLSIHVVTWNVASAAPPLDLSDLLQLNNRNLNLDIYVIGLQELNSGIIISLLSDAAFND SWSSFLMDVLSPLSFIKVSHVRMQGILLLVFAKYQHLPYIQLSTKSTPTGLFGYWGNKGGVNICLKLYG YYVSIINCHLPPHISNNYQRLEHFDRIEMQNCGRDIPNILDHDLIIFWGDMMNFRIEDFGLHVFRESIK NRCYGGLWEKDQLSIKKHDPDLLREFQEGRLFPPTYKFDRNSNDYDTSEKKRKPWTDRILWRLKRQPC AGPDTPIPPASHFSLSLRGYSSHMTYGISDHKPVSGTFDLELKPLVSAPLIVLMPEDLWTVENDMMVSY STSDFPSSPWDWIGLYKVGRLDVNDYVSYAWVGDSKVSCSDNLNQVYIDISNIPTTEDEFLCYSNSLR SVVGISRPFQIPPGSLREDPLGEAQPQI</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	50.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_057616](#)

Locus ID: 51763

UniProt ID: [Q9BT40](#)

RefSeq Size: 3001

Cytogenetics: 17p13.3

RefSeq ORF: 1344


Synonyms: MDCCAIID; PPS; SKIP

Summary: This gene encodes a protein with 5-phosphatase activity toward polyphosphate inositol. The protein localizes to the cytosol in regions lacking actin stress fibers. It is thought that this protein may negatively regulate the actin cytoskeleton. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Oct 2008]

Protein Families: Druggable Genome, Phosphatase

Protein Pathways: Inositol phosphate metabolism, Insulin signaling pathway, Metabolic pathways, Phosphatidylinositol signaling system

Product images:



Coomassie blue staining of purified INPP5K protein (Cat# TP300092). The protein was produced from HEK293T cells transfected with INPP5K cDNA clone (Cat# [RC200092]) using MegaTran 2.0 (Cat# [TT210002]).