

Product datasheet for **TP300095L**

CXorf26 (PBDC1) (NM_016500) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome X open reading frame 26 (CXorf26), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200095 protein sequence Red =Cloning site Green =Tags(s)
	<p>MAATSGTDEPVSGELVSVAHALSLPAESYGNDPDIEMAWAMRAMQHAEVYKLISSVDPQFLKLTQVDDQ IYSEFRKNFETLRIDVLDPEELKSESAKEKWRPFCLKFNGIVEDFNYGTLRLDSCQGYTEENTIFAPRI QFFAIEIARNREGYNKAVYISVQDKEGEGKGVNNGGKRAKSDGEEENTKNGGEGKADSGEEKEEGINREDK TDKGGEKGGKADKEINKSGEKAM</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	25.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_057584
Locus ID:	51260
UniProt ID:	Q9BVG4



[View online »](#)

RefSeq Size:	1214
Cytogenetics:	Xq13.3
RefSeq ORF:	699
Synonyms:	CXorf26

Product images:

Coomassie blue staining of purified PBDC1 protein (Cat# [TP300095]). The protein was produced from HEK293T cells transfected with PBDC1 cDNA clone (Cat# [RC200095]) using MegaTran 2.0 (Cat# [TT210002]).