

## Product datasheet for TP300130

### DNAJC17 (NM\_018163) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human Dnaj (Hsp40) homolog, subfamily C, member 17 (DNAJC17), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200130 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MAVTKELLQMDLYALLGIEEKAADKEVKKAYRQKALSCHPDKNPDNPRAAELFHQLSQALEVLTDAARA  
AYDKVRKAKKQAAERTQKLDEKRKKVKLDLEARERQAQAQSEEEEEESRSTRTLEQEIERLREEGSRQLE  
EQQLRIREQIRQERDQRLRGKAENTEGQGTPKLLKWKCKEDESKEGGYSKDVLLRLLQKYGEVNLNVLV  
SKKPGTAVVEFATVKAAELAVQNEVGLVDNPLKISWLEGQPQDAVGRSHSGLSKGSVLSERDYESLVMMR  
MRQAAERQQLIARMQEDQEGPPT

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	34.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_060633</a>
Locus ID:	55192



[View online »](#)

UniProt ID:	<u>Q9NVM6</u>
RefSeq Size:	1031
Cytogenetics:	15q15.1
RefSeq ORF:	912
Summary:	May negatively affect PAX8-induced thyroglobulin/TG transcription.[UniProtKB/Swiss-Prot Function]

### Product images:



Coomassie blue staining of purified DNAJC17 protein (Cat# TP300130). The protein was produced from HEK293T cells transfected with DNAJC17 cDNA clone (Cat# [RC200130]) using MegaTran 2.0 (Cat# [TT210002]).