

Product datasheet for **TP300163**

ZNF692 (NM_017865) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human zinc finger protein 692 (ZNF692), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200163 protein sequence Red =Cloning site Green =Tags(s)
	MASSPAVDVSCRRREKRRQLDARRSKCIRLGGHMEQWCLLKERLGFSLHSQLAKFLLDRYTSSGCVLCA GPEPLPPKGLQYLVLLSHAHSRECSLVPGLRGPGGQDGGVLWECSAGHTFSWGPSLSPTPSEAPKPASLP HTRRSWCSEATSGQELADLESEHDERTQEARLPRRVGPPPETFPPPGEIEGEEEDNDEDEEEMLSDAS LWTYSSPDSEPDAPRLLPSPVTCTPKEGETPPAPAALSSPLAVPALSASSLSSRAPPPAEVRVQPQLS RTPQAAQQTEALASTGSQAQSAPTPAWDEDTAQIGPKRIRKAAKRELMPCDFPGCGRIFSNRQYLNHHKK YQHIIHQKSFSCPEPACGKSFNFKKHLKEHMKLHSDTRDYICEFCARSFRTSSNLVIHRRHTGEKPLQCE ICGFTCRQKASLNWHQRKHAETVAALRFPCEFCGRFEKPDVAHAHRKSHPALLLAPQESPSGPLEPCP SISAPGPLGSSEGSRPSASPQAPTLLPQQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	56.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq:	<u>NP_060335</u>
Locus ID:	55657
UniProt ID:	<u>Q9BU19</u>
RefSeq Size:	1992
Cytogenetics:	1q44
RefSeq ORF:	1557
Synonyms:	AREBP; Zfp692
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified ZNF692 protein (Cat# TP300163). The protein was produced from HEK293T cells transfected with ZNF692 cDNA clone (Cat# [RC200163]) using MegaTran 2.0 (Cat# [TT210002]).