

## Product datasheet for TP300170M

### VSIG2 (NM\_014312) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human V-set and immunoglobulin domain containing 2 (VSIG2), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200170 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MAELPGPFLCGALLGFLCLSGLAVEVKVPTEPLSTPLGKTAELTCTYSTSVGDSFALEWSFVQPGKPISE  
SHPILYFTNGHLYPTGSKSKRVSLQNPPTVGVATLKLTDVHPSDTGTYLQVNNPPDFYTNGLGLINLT  
VLVPPSNPLCSQSQGQTSVGGSTALRCSSEGAPKPVYNWVRLGTFPTSPGSMVQDEVSGQLILTNSLT  
SSGTYRCVATNQMGASCELTLVTEPSQGRVAGALIGVLLGVLLLSVAAFCLVRFQKERGKKPKETYGG  
SDLREDAIAPGISEHTCMRADSSKGFLERPSSASTVTTTTSKSLPMVV

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**


Tag:	C-Myc/DDK
Predicted MW:	34.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_055127</a>
Locus ID:	23584



[View online »](#)

UniProt ID:	<a href="#">Q96IQ7</a>
RefSeq Size:	1138
Cytogenetics:	11q24.2
RefSeq ORF:	981
Synonyms:	2210413P10Rik; CTH; CTXL
Protein Families:	Druggable Genome, Transmembrane

### Product images:



Coomassie blue staining of purified VSIG2 protein (Cat# [TP300170]). The protein was produced from HEK293T cells transfected with VSIG2 cDNA clone (Cat# [RC200170]) using MegaTran 2.0 (Cat# [TT210002]).