

Product datasheet for TP300175L

DUS2L (DUS2) (NM_017803) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human dihydrouridine synthase 2-like, SMM1 homolog (<i>S. cerevisiae</i>) (DUS2L), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200175 protein sequence Red=Cloning site Green=Tags(s)

MILNSLSLCYHNKILAPMVRVGTLP MRLALDYGADIVYCEELIDLKMIQCKRVNEVLSTVDFVAPDD
RVVFRTCEREQNRVVFQMGTSDAERALAVARLVENDVAGIDVNMGC PKQYSTKGGMGAALLSDPKIEKI
LSTLVKGTTRPV TCKIRILPSLED TSLVKRIERTGIAAIVHGRKREERPQHPVSCEVIKAIAD TLSIP
VIANGGSHDHIQQYSDIEDFRQATAASSVMVARAAMWNPSIFLKEGLRPLEEVMQKYIRYAVQYDNHYTN
TKYCLCQMLREQL ESPQGRLLHAAQSSREICEAFGLGAFYEETTQELDAQQARLSAKTSEQTGEP AEDTS
GVIKMAVKFDRRAYPAQITPKMCLLEWCRREKLAQPVYETVQRPLDRL FSSIVTVAEQKYQSTLWDKSKK
LAEQAAAIVCLRSQGLPEGRLGEESPSLHKRKREAPDQDPGGPRAQELAQPGDLCKKPFVALGSGEESPL
EGW

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	54.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_060273
Locus ID:	54920
UniProt ID:	Q9NX74
RefSeq Size:	2342
Cytogenetics:	16q22.1
RefSeq ORF:	1479
Synonyms:	DUS2L; SMM1; URLC8
Summary:	This gene encodes a cytoplasmic protein that catalyzes the conversion of uridine residues to dihydrouridine in the D-loop of tRNA. The resulting modified bases confer enhanced regional flexibility to tRNA. The encoded protein may increase the rate of translation by inhibiting an interferon-induced protein kinase. This gene has been implicated in pulmonary carcinogenesis. Alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Nov 2012]

Product images:



Coomassie blue staining of purified DUS2 protein (Cat# [TP300175]). The protein was produced from HEK293T cells transfected with DUS2 cDNA clone (Cat# [RC200175]) using MegaTran 2.0 (Cat# [TT210002]).