


Locus ID:	55603
UniProt ID:	Q96IP4
RefSeq Size:	5617
Cytogenetics:	6q14.1
RefSeq ORF:	1326
Synonyms:	C6orf37; FAM46A; OI18; XTP11
Summary:	Probable nucleotidyltransferase that may act as a non-canonical poly(A) RNA polymerase. [UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified TENT5A protein (Cat# TP300181). The protein was produced from HEK293T cells transfected with TENT5A cDNA clone (Cat# [RC200181]) using MegaTran 2.0 (Cat# [TT210002]).