

## Product datasheet for TP300208M

### MCM3 (NM\_002388) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human minichromosome maintenance complex component 3 (MCM3), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200208 protein sequence Red=Cloning site Green=Tags(s)

MAGTVLDDVELREAQRDYLDLDFLDEEDQGIYQSKVRELISDNQYRLIVNVNDLRRKNEKRANLLNNAF  
EELVAFQRALKDFVASIDATYAKQYEEFYVGLGSGFKHVSPTLTSCFLSCVVCVEGIVTKCSLVRPK  
VRSVHYCPATKKTERRYSIDLTLVAFPSVYPTKDEENPLETEYGLSVYKDHQTITIQEMPEKAPA  
GQLPRSDVILDDDLVDKAKPGDRVQVVGTYRCLPGKGGYTSCTFRFTVLIACNVKQMSKDAQPSFSAED  
IAKIKKFSKTRSKDIFDQLAKSLAPSIHGHDYVKKAILCLLLGGVERDLENGSHIRGDINILLIGDPSVA  
KSQLLRYVLTAPRAIPTTGRGSSGVGLTAAVTTDQETGERRLEAGAMVLADRGVVCIDEDFKMSDMDRT  
AIHEVMEQGRVTIAKAGIHARLNARCSVLAANPVYGRYDQYKTPMENIGLQDSSLRFDLLFIMLDQMD  
PEQDREISDHVLRMHRYRAPGEQDGDAMPLGSAVDILATDDPNFSQEDQQDQTIYEKHDNLLHGTTKKKE  
KMVSAAFMKKYIHVAKIIPVLTQESATYIAEYSRLRSQDSMSSDTARTSPVTARTLETIIRLATAHAK  
ARMSKTVDLQDAEEAVELVQYAYFKKVEKEKKRKRSEDESETEDEEEKSQEDQEQRKRKRTRQPDAK  
DGDSYDPYDFSDTEEMPQVHTPKTADSQETKESQKVELSESRLKAFKVALLDVFREAHAQSIGMNRLTE  
SINRDSEEPFSSVEIQAALSKMQDDNQVMVSEGIIFLI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	90.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

**Note:** For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

**Storage:** Store at -80°C.

**Stability:** Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

**RefSeq:** [NP\\_002379](#)

**Locus ID:** 4172

**UniProt ID:** [P25205](#), [A0A0S2Z4T1](#)

**RefSeq Size:** 3234

**Cytogenetics:** 6p12.2

**RefSeq ORF:** 2424


**Synonyms:** HCC5; P1-MCM3; P1.h; RLFB

**Summary:** The protein encoded by this gene is one of the highly conserved mini-chromosome maintenance proteins (MCM) that are involved in the initiation of eukaryotic genome replication. The hexameric protein complex formed by MCM proteins is a key component of the pre-replication complex (pre\_RC) and may be involved in the formation of replication forks and in the recruitment of other DNA replication related proteins. This protein is a subunit of the protein complex that consists of MCM2-7. It has been shown to interact directly with MCM5/CDC46. This protein also interacts with and is acetylated by MCM3AP, a chromatin-associated acetyltransferase. The acetylation of this protein inhibits the initiation of DNA replication and cell cycle progression. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2018]

**Protein Families:** Druggable Genome, Stem cell - Pluripotency, Transcription Factors

**Protein Pathways:** Cell cycle, DNA replication

### Product images:



Coomassie blue staining of purified MCM3 protein (Cat# [TP300208]). The protein was produced from HEK293T cells transfected with MCM3 cDNA clone (Cat# [RC200208]) using MegaTran 2.0 (Cat# [TT210002]).