

Product datasheet for **TP300264L**

MSF (SEPT9) (NM_006640) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human septin 9 (SEPT9), transcript variant 3, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200264 protein sequence Red =Cloning site Green =Tags(s)
	<p>MERDRISALKRSFEVEEVETPNSTPPRRVQTPLLRTVASSTQKFQDLGVKNSEPSARHVDSLSQRSPKA SLRRVELSGPKAAEPVSRTELSIDISSKQVENAGAIGPSRFLKRAEVLGHKTPEPAPRRTEITIVKPQ ESAHRREPPASKVPEVPTAPATDAAPKRVEIQMPKPAEAPTAPSPAQTLENSEPAPVSQLQSRLEPKPQ PPVAEATPRSQEATEAAPSCVGDMDTPRDAGLKQAPASRNEKAPVDFGYVGIDSIQEMRRKAMKQGFE FNIMVVGQSGLGKSTLINTLFKSKSRKSVQPTSEERIPKTIEIKSITHDIEEKGVRMKLTVIDTPGFGD HINNENCWQPIMKFINDQYEKYLQEEVNINRKKRIPDTRVHCCLYFIPATGHSLRPLDIEFMKRLSKVNV IVPIAKADTLTLEERVHFKQRITADLLSNGIDVYPQKEFDEDESDRLVNEKFREMIPFAVVGSDHEYQV NGKRILGRKTKWGTIEVENTTHCEFAYLRDLLIRTHMQNIKDITSSIHFEAYRVKRLNEGSSAMANGVEE KEPEAPEM</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	63.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_006631
Locus ID:	10801
UniProt ID:	Q9UHD8 , A0A0S2Z5A5
RefSeq Size:	4469
Cytogenetics:	17q25.3
RefSeq ORF:	1704
Synonyms:	AF17q25; MSF; MSF1; NAPB; PNUTL4; SEPT9; SeptD1; SINT1
Summary:	This gene is a member of the septin family involved in cytokinesis and cell cycle control. This gene is a candidate for the ovarian tumor suppressor gene. Mutations in this gene cause hereditary neuralgic amyotrophy, also known as neuritis with brachial predilection. A chromosomal translocation involving this gene on chromosome 17 and the MLL gene on chromosome 11 results in acute myelomonocytic leukemia. Multiple alternatively spliced transcript variants encoding different isoforms have been described.[provided by RefSeq, Mar 2009]
Protein Families:	Druggable Genome

Product images:


Coomassie blue staining of purified SEPTIN9 protein (Cat# [TP300264]). The protein was produced from HEK293T cells transfected with SEPTIN9 cDNA clone (Cat# [RC200264]) using MegaTran 2.0 (Cat# [TT210002]).