

Product datasheet for **TP300267**

COPS8 (NM_198189) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human COP9 constitutive photomorphogenic homolog subunit 8 (Arabidopsis) (COPS8), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200267 protein sequence Red =Cloning site Green =Tags(s)
	 MPVAVMAESAFSFKLLDQCENQELEAPGGIATPPVYGQLLALYLLHNDMNNARYLWKRIIPPAIKSANSE LGGIWSVGQRIWQRDFPGIYTTINAHQWSETVQPIMEALRDATRRRAFALVSQAYTSIIADDFAAVGLP VEEAVKGILEQGWQADSTTRMVLPRKPVAGALDVSFNKFIPLSEPAVPPPIPNEQQLARLTDYVAFLEN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	17.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_937832
Locus ID:	10920
UniProt ID:	Q99627



[View online »](#)

RefSeq Size: 2339


Cytogenetics: 2q37.3

RefSeq ORF: 630

Synonyms: COP9; CSN8; SGN8

Summary: The protein encoded by this gene is one of the eight subunits of COP9 signalosome, a highly conserved protein complex that functions as an important regulator in multiple signaling pathways. The structure and function of COP9 signalosome is similar to that of the 19S regulatory particle of 26S proteasome. COP9 signalosome has been shown to interact with SCF-type E3 ubiquitin ligases and act as a positive regulator of E3 ubiquitin ligases. Alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified COPS8 protein (Cat# TP300267). The protein was produced from HEK293T cells transfected with COPS8 cDNA clone (Cat# [RC200267]) using MegaTran 2.0 (Cat# [TT210002]).