

Product datasheet for **TP300288L**

LSM1 (NM_014462) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human LSM1 homolog, U6 small nuclear RNA associated (S. cerevisiae) (LSM1), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200288 protein sequence Red =Cloning site Green =Tags(s)
	MNYMPGTASLIEDIDKKHLVLLRDGRTLIGFLRSIDQFANLVLHQTVERIHVGKKYGDIPRGIFVVRGEN VLLGEIDLEKESDTPLQQVSIEEILEEQRVEQQTKLEAEKLVQALKDRGLSIPRADTLDEY
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	15 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_055277</u>
Locus ID:	27257
UniProt ID:	<u>O15116</u> , <u>A0A0S2Z590</u>
RefSeq Size:	1161



[View online »](#)

Cytogenetics: 8p11.23

RefSeq ORF: 399


Synonyms: CASM; YJL124C

Summary: This gene encodes a member of the LSM family of RNA-binding proteins. LSM proteins form stable heteromers that bind specifically to the 3'-terminal oligo(U) tract of U6 snRNA and may play a role in pre-mRNA splicing by mediating U4/U6 snRNP formation. Increased expression of this gene may play a role in cellular transformation and the progression of several malignancies including lung cancer, mesothelioma and breast cancer. Alternatively spliced transcript variants have been observed for this gene, and a pseudogene of this gene is located on the short arm of chromosome 9. [provided by RefSeq, Nov 2011]

Protein Families: Stem cell - Pluripotency

Protein Pathways: RNA degradation

Product images:



Coomassie blue staining of purified LSM1 protein (Cat# [TP300288]). The protein was produced from HEK293T cells transfected with LSM1 cDNA clone (Cat# [RC200288]) using MegaTran 2.0 (Cat# [TT210002]).