

Product datasheet for **TP300290L**

TFIP11 (NM_001008697) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human tuftelin interacting protein 11 (TFIP11), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200290 protein sequence Red =Cloning site Green =Tags(s)

MSLSHLYRDGEGRIDDDDDERENFEITDWDLQNEFNPNRQRHWQTKKEATYGVWAERDSDDERPSFGGKR
ARDYSAPVNFISAGLKKGAEEAELESDDEEKPVKQDDFPKDFGPRKLKTGGNFKPSQKGFAGGTKSFM
DFGSWERHTKGIGQKLLQKMGYVPRGLGKNAQGIINPIEAKQRKGGAVGAYGSERTTQSMQDFPVVDS
EEEEEEEFQKELSQRKDPGSGSKKPKYSYKVEELKAKGRISKKLTAPQKELSQVKVIDMTGREQKVYY
SYSQISHKHNPDDGLPLQSQQLPQSGKEAKAPGFALPELEHNLQLLIDLTEQEIIQNDRQLQYERDMVV
NLFHELEKMTEVLDHEERVISNLSKVLEMVEECERRMQPDCSNPLTLDECARIFETLQDKYEEYRMSDR
VDLAVAIVYPLMKEYFKWDPLKDCYGTETIISKWKSLENDQLLSHGGQDLSADAFHRLIWEVWMPFVR
NIVTQWQPRNCDPMVDFLDSWVHIIPWILDNILDQLIFPKLQKEVENWNPLTDTVPIHSHWHPWLPLMQ
ARLEPLYSPIRSKLSSALQKWHPSDSSAKLILQPWKDVFTPGSWEAFMVKNIVPKLGMCLGELVINPHQQ
HMDAFYWVIDWEGMISVSSLVGLLEKHFFPKWLQVLCWSLSNSPNYEEITKWYLGWKSMSFSDQVLAHPSV
KDKFNEALDIMNRAVSSNVGAYMQPGARENIAYLHTERRKDFQYEAMQERREAENMAQRGIGVAASSVP
MNFKDLIETKAEHNIVFMPVIGKRHEGKQLYTFGRIVYIDRGVVFVQGEKTWVPTSLQSLIDMAK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	96.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001008697
Locus ID:	24144
UniProt ID:	Q9UBB9 , A0A024R117
RefSeq Size:	2941
Cytogenetics:	22q12.1
RefSeq ORF:	2511
Synonyms:	bK445C9.6; NTR1; Spp382; STIP; STIP-1; TIP39
Summary:	This gene encodes a protein component of the spliceosome that promotes the release of the lariat-intron during late-stage splicing through the recruitment of a pre-mRNA splicing factor called DEAH-box helicase 15. The encoded protein contains a G-patch domain, a hallmark of RNA-processing proteins, that binds DEAH-box helicase 15. This protein contains an atypical nuclear localization sequence as well as a nuclear speckle-targeting sequence, enabling it to localize to distinct speckled regions within the cell nucleus. Polymorphisms in this gene are associated with dental caries suggesting a role in amelogenesis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2016]

Product images:



Coomassie blue staining of purified TFIP11 protein (Cat# [TP300290]). The protein was produced from HEK293T cells transfected with TFIP11 cDNA clone (Cat# [RC200290]) using MegaTran 2.0 (Cat# [TT210002]).