

Product datasheet for **TP300300M**

GPIP137 (CAPRIN1) (NM_005898) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human cell cycle associated protein 1 (CAPRIN1), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200300 protein sequence Red =Cloning site Green =Tags(s) MPSATSHSGSGSKSSGPPPPSGSSGSEAAAGAGAAAPASQHPATGTGAVQTEAMKQILGVIDKKLRNLEK KKGKLDYQERMNKGRLNQDQLDAVSKYQEVNTNNLEFAKELQRSFMALSQDIQTIKKTARREQLMREE AEQKRLKTVLELQYVLDKLGDDDEVRTDLKQGLNGVPILSEEELSLLDDEFYKLVDPERDMSLRLENEQEHA SIHLWDLLEGKEKPVCGTTYKVLKEIVERVFQSNYFDSTHNHQNGLCSEEEEAASAPAVEDQVPEAEPEPA EEYTEQSEVESTYVNRQFMAETQFTSGEKEQVDEWTVETVEVNSLQQPQAASPSVPEPHSLTPVAQA DPLVRRQRVQDLMAQMQGPYNFIQDSMLDFENQTLDPAIVSAQPMNPTQNMMPQLVCPVHSESRLAQP NQVPVQPEATQVPLVSSTSEGYTASQPLYQPSHATEQRPOKEPIDQIATISLNTDQTTASSSLPAASQP QVFQAGTSKPLHSSGINVNAAPFQSMQTVFNMNAPVPPVNEPETLKQQNQYQASYNQSFSSQPHQVEQTE LQQEQLQTVVGTYHGSPDQSHQVTGNHQPPQNTGFPRSNQPYNSRGVSRGGSRGARGLMNGYRGPAN GFRGGYDGYRPSFSNTPNSGYTQSQFSAPRDYSGYQRDGYQQNFKRGSQSGPRGAPRGRGGPPRPNRGM PQMNTQQVN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	78.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Bioactivity:	Co-immunoprecipitation (PMID: 25784705)
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_005889](#)

Locus ID: 4076

UniProt ID: [Q14444](#)

RefSeq Size: 5562


Cytogenetics: 11p13

RefSeq ORF: 2127


Synonyms: GPIAP1; GPIP137; GRIP137; M11S1; p137GPI; RNG105

Summary: May regulate the transport and translation of mRNAs of proteins involved in synaptic plasticity in neurons and cell proliferation and migration in multiple cell types. Binds directly and selectively to MYC and CCND2 RNAs. In neuronal cells, directly binds to several mRNAs associated with RNA granules, including BDNF, CAMK2A, CREB1, MAP2, NTRK2 mRNAs, as well as to GRIN1 and KPNB1 mRNAs, but not to rRNAs.[UniProtKB/Swiss-Prot Function]

Product images:



Caprin1 interacts with inactive PKR. Purified Caprin1 (OriGene [TP300300]) was incubated with inactive PKR (PKR) or active PKR (P-PKR), immunoprecipitated with either nonspecific IgG or caprin1 antibodies, and analyzed by Western blotting for caprin1 or total PKR with corresponding antibodies. Figure cited from MBio, PMID: 25784705



Coomassie blue staining of purified CAPRIN1 protein (Cat# [TP300300]). The protein was produced from HEK293T cells transfected with CAPRIN1 cDNA clone (Cat# [RC200300]) using MegaTran 2.0 (Cat# [TT210002]).