

## Product datasheet for **TP300303M**

### WDR1 (NM\_017491) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human WD repeat domain 1 (WDR1), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC200303 protein sequence
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MPYEIKKVFASLPQVERGVSKIIGGDPKGNFLYTNGKCVILRNIDNPALADIYTEHAHQVWAKYAPSG  
FYIASGDVSGKLRWDTTQKEHLLKYEYQPFAGKIKDIAWTEDESKRIAVVGEKREKFGAVFLWDSGSSVG  
EITGHNVINSVDIKQSRPYRLATGSDDNCAAFFEGPPFKFTIGDHSRVNVCVRFSPDGNRFATASAD  
GQIYIDGKTGEKVCALGGSKAHDGGIYAIWSPDSTHLLSASGDKTSKIWDVSVNSVSTFPMGSTVLD  
QQLGCLWQKDHLLSVLSGYINYLDNRNPSKPLHVIKGHSKSIQCLTVHKNGGKSYIYSGSHDGHINYWD  
SETGENDSFAGKGTNQVSRMTVDESGQLISCSMDDTVRYTSLMLRDYSGQGKLDVQPKCVAVGPGGY  
AVVVICIGQIVLLKDKRKCFSIDNPGYEPEVAVHPGGDTVAIGGVDGNVRLYSILGTTLKDEGKLLLEAKG  
PVTDVAYSHDGAFLAVCDASKVWTVFSVADGYSENNVFGHHAKIVCLAWSPDNEHFASGGMDMMVYVWT  
LSDPETRVKIQDAHRLHHVSSLAWLDEHTLVTTSHDASVKEWTITY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	66 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_059830</a>
<b>Locus ID:</b>	9948
<b>UniProt ID:</b>	<a href="#">O75083</a> , <a href="#">V9HWG7</a>
<b>RefSeq Size:</b>	3160
<b>Cytogenetics:</b>	4p16.1
<b>RefSeq ORF:</b>	1818
<b>Synonyms:</b>	AIP1; HEL-S-52; NORI-1; PFITS
<b>Summary:</b>	This gene encodes a protein containing 9 WD repeats. WD repeats are approximately 30- to 40-amino acid domains containing several conserved residues, mostly including a trp-asp at the C-terminal end. WD domains are involved in protein-protein interactions. The encoded protein may help induce the disassembly of actin filaments. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

### Product images:



Coomassie blue staining of purified WDR1 protein (Cat# [TP300303]). The protein was produced from HEK293T cells transfected with WDR1 cDNA clone (Cat# [RC200303]) using MegaTran 2.0 (Cat# [TT210002]).