

Product datasheet for **TP300343L**

UXT (NM_004182) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ubiquitously-expressed transcript (UXT), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200343 protein sequence Red =Cloning site Green =Tags(s)

MATPPKRRAVEATGEKVLRYETFISDVLQRDLRKVLDHRDKVYEQLAKYLQLRNVIERLQEAKHSELYMQ
VDLGCNFFVDTVVPDTSRIYVALGYGFFLELTLAEALKFIDRKSSLLTELSNSLTKDSMNIKAHIHMLLE
GLRELQGLQNFPEKPHH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	18.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_004173</u>
Locus ID:	8409
UniProt ID:	<u>Q9UBK9</u>



[View online »](#)

RefSeq Size: 619

Cytogenetics: Xp11.23


RefSeq ORF: 471

Synonyms: ART-27; STAP1

Summary: The protein encoded by this gene functions as a cofactor that modulates androgen receptor-dependent transcription, and also plays a critical role in tumor necrosis factor-induced apoptosis. Expression of this gene may play a role in tumorigenesis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Nov 2011]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified UXT protein (Cat# [TP300343]). The protein was produced from HEK293T cells transfected with UXT cDNA clone (Cat# [RC200343]) using MegaTran 2.0 (Cat# [TT210002]).