

## Product datasheet for **TP300449M**

### Nck beta (NCK2) (NM\_001004720) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human NCK adaptor protein 2 (NCK2), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200449 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MTEEVIVIAKWDYTAQQDQELDIKKNERLWLLDDSKTWWRVNRNAANRTGYVPSNYVERKNSLKKGSLVKN LKDTLGLGKTRRKSARDASPTPSTDAEYPANGSGADRIYDLNIPAFVKFAYVAEREDELVLKGSRVTV MEKCSGWWRGSYNGQIGWFPSNYVLEEVDEAAAESPSFLSLRKASLSNGQGSRVLHVVTLYPFSSVT EEELNFEKGETMEVIEKPENDEPWWKCKNARGQVGLVPKNYVWLSDGPALHPAHAPQISYTGPSSSGRF AGREWYYGNVTRHQAECALNERGVEGDFLIRDSESSPSDFSVSLKASGKNKHFVKQLVDNVYICIGQRRFH TMDELVEHYKKAPIFTSEHGKLYLVRALQ
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	42.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_001004720</a></u>
Locus ID:	8440



[View online »](#)

UniProt ID: [Q43639](#), [A0A0S2Z4M6](#)

RefSeq Size: 2417

Cytogenetics: 2q12.2

RefSeq ORF: 1140


Synonyms: GRB4; NCKbeta

**Summary:** This gene encodes a member of the NCK family of adaptor proteins. The protein contains three SH3 domains and one SH2 domain. The protein has no known catalytic function but has been shown to bind and recruit various proteins involved in the regulation of receptor protein tyrosine kinases. It is through these regulatory activities that this protein is believed to be involved in cytoskeletal reorganization. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]

**Protein Families:** Druggable Genome

**Protein Pathways:** Axon guidance, ErbB signaling pathway, Pathogenic Escherichia coli infection, T cell receptor signaling pathway

### Product images:



Coomassie blue staining of purified NCK2 protein (Cat# [TP300449]). The protein was produced from HEK293T cells transfected with NCK2 cDNA clone (Cat# [RC200449]) using MegaTran 2.0 (Cat# [TT210002]).