

Product datasheet for **TP300468L**

GSK3 beta (GSK3B) (NM_002093) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human glycogen synthase kinase 3 beta (GSK3B), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200468 representing NM_002093 Red =Cloning site Green =Tags(s)
	 MSGRPRTTSAESCKPVQQPSAFGSMKVS RDKDGSKVTTVATPGQGPDRPQEVSYTDTKVI GNNGSFGW YQAKLCDSGELVAIKKVLQDKRFKNRELQIMRKLDHCNIVRLRYFFYSSGEKKDEVYLNVLVDYVPETVY RVARHYSRAKQTLPIYVKLYMYQLFRSLAYIHSFGICH RDIKPQNLLLDPDTAVLKLCDFGSAKQLVRG EPNVSYICSRYYRAPELIFGATDYTSSIDVWSAGCVLAELLGQPIFGDSGVDQLVEIKVLGTPREQ IREMNPNYTEFKFPQIKAHPWTKDSSGTGHFTSGVRVFRPRTPPEAIALCSRLLEYTP TARLTPLEACAH SFFDEL RDPNVKLPNGRDTPALFNFTTQELSSNPPLATILIPPHARIQAAASTPTNATAASDANTGDRGQ TNNAASASASNST TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	47.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_002084</u>



[View online »](#)


Locus ID: 2932
UniProt ID: [P49841](#)
RefSeq Size: 1639
Cytogenetics: 3q13.33
RefSeq ORF: 1299

Summary: The protein encoded by this gene is a serine-threonine kinase belonging to the glycogen synthase kinase subfamily. It is a negative regulator of glucose homeostasis and is involved in energy metabolism, inflammation, ER-stress, mitochondrial dysfunction, and apoptotic pathways. Defects in this gene have been associated with Parkinson disease and Alzheimer disease. [provided by RefSeq, Aug 2017]

Protein Families: Druggable Genome, Protein Kinase

Protein Pathways: Alzheimer's disease, Axon guidance, Basal cell carcinoma, B cell receptor signaling pathway, Cell cycle, Chemokine signaling pathway, Colorectal cancer, Endometrial cancer, ErbB signaling pathway, Focal adhesion, Hedgehog signaling pathway, Insulin signaling pathway, Melanogenesis, Neurotrophin signaling pathway, Pathways in cancer, Prostate cancer, T cell receptor signaling pathway, Wnt signaling pathway

Product images:



Coomassie blue staining of purified GSK3B protein (Cat# [TP300468]). The protein was produced from HEK293T cells transfected with GSK3B cDNA clone (Cat# [RC200468]) using MegaTran 2.0 (Cat# [TT210002]).