

## Product datasheet for TP300519L

### GFPT2 (NM\_005110) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human glutamine-fructose-6-phosphate transaminase 2 (GFPT2), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200519 protein sequence Red=Cloning site Green=Tags(s)

MCGIFAYMNYRVPRTTRKEIFETLIKGLQRLEYRGYDSAGVAIDGNNHEVKERHIQLVKKRGKVKALDEEL  
YKQDSMDLKVEFETHFGIAHTRWATHGVPSAVNSHPQRSKGNFVVIHNGIITNYKDLRKFLESKGYEF  
ESETDTETIAKLIKVFVDNRETEDITFSTLVERVIQQLLEGAFALVFKSVHYPGEAVATRRGSPLLIGVRS  
KYKLSTEQIPILYRTCTLENVKNICKTRMKRLDSSACLHAVGDKAVEFFASDASAIIEHTNRVIFLEDD  
DIAAVADGKLSIHRVKRSASDDPSRAIQTLMELQQIMKGNFSAFMQKEIFEQPESVFNTMRGRVNFETN  
TVLLGGLKDHLKEIRRCRLVIGCGTSYHAAVATRQVLEELTELPVMVELASDFLDRNTPVFRDDVCFE  
ISQSGETADTLLALRYCKDRGALTGVGTNTVGVSSISRETDCGVHINAGPEVGVASTKAYTSQFISLVMFG  
LMMSEDRISLQNRQEIIRGLRSLPELIKEVLSLEEKIHDLALEYTQRSLVMGRGYNATCLEGALKI  
KEITYMHSEGILAGELKHGPLALIDKQMPVIMVIMKDPCKAFKQNALQQVTARQGRPIILCSKDDTTESSK  
FAYKTIELPHTVDCLQGILSVIPLQLLSFHLAVLRGYDVFPRNLAKSVTVE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	76.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_005101</a>
<b>Locus ID:</b>	9945
<b>UniProt ID:</b>	<a href="#">O94808</a> , <a href="#">A0A0S2Z4X9</a> , <a href="#">B3KMR8</a>
<b>RefSeq Size:</b>	3011
<b>Cytogenetics:</b>	5q35.3
<b>RefSeq ORF:</b>	2046
<b>Synonyms:</b>	GFAT; GFAT 2; GFAT2
<b>Summary:</b>	Controls the flux of glucose into the hexosamine pathway. Most likely involved in regulating the availability of precursors for N- and O-linked glycosylation of proteins.[UniProtKB/Swiss-Prot Function]
<b>Protein Families:</b>	Protease
<b>Protein Pathways:</b>	Alanine, aspartate and glutamate metabolism, Amino sugar and nucleotide sugar metabolism, Metabolic pathways

### Product images:



Coomassie blue staining of purified GFPT2 protein (Cat# [TP300519]). The protein was produced from HEK293T cells transfected with GFPT2 cDNA clone (Cat# [RC200519]) using MegaTran 2.0 (Cat# [TT210002]).