

## **Product datasheet for TP300559**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **HOXA9 (NM\_152739) Human Recombinant Protein**

**Product data:** 

**Product Type:** Recombinant Proteins

**Description:** Recombinant protein of human homeobox A9 (HOXA9), 20 μg

Species: Human
Expression Host: HEK293T

**Expression cDNA Clone** >RC200559 representing NM\_152739 **or AA Sequence:** Red=Cloning site Green=Tags(s)

MATTGALGNYYVDSFLLGADAADELSVGRYAPGTLGQPPRQAATLAEHPDFSPCSFQSKATVFGASWNPV HAAGANAVPAAVYHHHHHHPYVHPQAPVAAAAPDGRYMRSWLEPTPGALSFAGLPSSRPYGIKPEPLSAR RGDCPTLDTHTLSLTDYACGSPPVDREKQPSEGAFSENNAENESGGDKPPIDPNNPAANWLHARSTRKKR

CPYTKHQTLELEKEFLFNMYLTRDRRYEVARLLNLTERQVKIWFQNRRMKMKKINKDRAKDE

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV** 

Tag: C-Myc/DDK

Predicted MW: 30 kDa

**Concentration:** >0.05 μg/μL as determined by microplate BCA method

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining

**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

**Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by conventional

chromatography steps.

**Note:** For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

**RefSeq:** NP 689952

**Locus ID:** 3205

UniProt ID: P31269





RefSeq Size: 2076

Cytogenetics: 7p15.2 RefSeq ORF: 816


**Synonyms:** ABD-B; HOX1; HOX1.7; HOX1G

**Summary:** In vertebrates, the genes encoding the class of transcription factors called homeobox genes

are found in clusters named A, B, C, and D on four separate chromosomes. Expression of these proteins is spatially and temporally regulated during embryonic development. This gene is part of the A cluster on chromosome 7 and encodes a DNA-binding transcription factor which may regulate gene expression, morphogenesis, and differentiation. This gene is highly similar to the abdominal-B (Abd-B) gene of Drosophila. A specific translocation event which causes a fusion between this gene and the NUP98 gene has been associated with myeloid leukemogenesis. Read-through transcription exists between this gene and the upstream

homeobox A10 (HOXA10) gene.[provided by RefSeq, Mar 2011]

## **Product images:**



Coomassie blue staining of purified HOXA9 protein (Cat# TP300559). The protein was produced from HEK293T cells transfected with HOXA9 cDNA clone (Cat# [RC200559]) using MegaTran 2.0 (Cat# [TT210002]).