

## **Product datasheet for TP300583M**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## ID3 (NM\_002167) Human Recombinant Protein

**Product data:** 

**Product Type:** Recombinant Proteins

**Description:** Recombinant protein of human inhibitor of DNA binding 3, dominant negative helix-loop-helix

protein (ID3), 100 µg

Species: Human
Expression Host: HEK293T

**Expression cDNA Clone** >RC200583 protein sequence or AA Sequence: Red=Cloning site Green=Tags(s)

MKALSPVRGCYEAVCCLSERSLAIARGRGKGPAAEEPLSLLDDMNHCYSRLRELVPGVPRGTQLSQVEIL

QRVIDYILDLQVVLAEPAPGPPDGPHLPIQTAELAPELVISNDKRSFCH

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV** 

Tag: C-Myc/DDK

**Predicted MW:** 12.8 kDa

Concentration:  $>0.05 \mu g/\mu L$  as determined by microplate BCA method

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining

**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

**Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by

conventional chromatography steps.

**Note:** For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 002158

 Locus ID:
 3399

 UniProt ID:
 Q02535

 RefSeq Size:
 1252





Cytogenetics: 1p36.12

RefSeq ORF: 357

Synonyms: bHLHb25; HEIR-1

**Summary:** The protein encoded by this gene is a helix-loop-helix (HLH) protein that can form


> heterodimers with other HLH proteins. However, the encoded protein lacks a basic DNAbinding domain and therefore inhibits the DNA binding of any HLH protein with which it

interacts. [provided by RefSeq, Aug 2011]

**Protein Families:** Druggable Genome, ES Cell Differentiation/IPS, Transcription Factors

**Protein Pathways:** TGF-beta signaling pathway

## **Product images:**



Coomassie blue staining of purified ID3 protein (Cat# [TP300583]). The protein was produced from HEK293T cells transfected with ID3 cDNA clone (Cat# [RC200583]) using MegaTran 2.0

(Cat# [TT210002]).