

Product datasheet for **TP300589L**

Calumenin (CALU) (NM_001219) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human calumenin (CALU), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200589 protein sequence Red =Cloning site Green =Tags(s)
	 MDLRQFLMCLSLCTAFALSKPTEKKDRVHHEPQLSDKVHNDAAQSFYDHDHDAFLGAEAAKTFDQLTPESK ERLGKIVSKIDGDKDGFVTVDELKDWIKFAQKRWIYEDVERQWKGHDLNEDGLVSWEEYKNATYGYVLDL PDPDDGFNYKQMMVRDERRFKMADKDGDLIATKEEFTAFLHPEEYDYMKDIVVQETMEDIDKNADGFIDL EEYIGDMYSHDGNTDEPEWVKTEREQFVEFRDKNRDGKMDKEETKDWILPSDYDHAEAEARHLVYESDQN KDGLTKKEIVDKYDLFVGSQATDFGEALVRHDEF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	34.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001210</u>
Locus ID:	813



[View online »](#)

UniProt ID: [O43852](#), [Q6IAW5](#), [B3KQF5](#), [D6QS48](#)

RefSeq Size: 5366


Cytogenetics: 7q32.1

RefSeq ORF: 945

Summary: The product of this gene is a calcium-binding protein localized in the endoplasmic reticulum (ER) and it is involved in such ER functions as protein folding and sorting. This protein belongs to a family of multiple EF-hand proteins (CERC) that include reticulocalbin, ERC-55, and Cab45 and the product of this gene. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Oct 2008]

Protein Families: Secreted Protein

Product images:



Coomassie blue staining of purified CALU protein (Cat# [TP300589]). The protein was produced from HEK293T cells transfected with CALU cDNA clone (Cat# [RC200589]) using MegaTran 2.0 (Cat# [TT210002]).