

Product datasheet for **TP300594M**

NMT1 (NM_021079) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human N-myristoyltransferase 1 (NMT1), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200594 representing NM_021079 Red =Cloning site Green =Tags(s)
	<p>MADESETAVKPPAPPLPQMMEGNGNGHEHCSDCENEEDNSYNRGGSPANDTGAKKKKKKQKKKKEKGSE TDSAQDQPVKMNSLPAERIQEIQKAIELFSVGQGPACTMEEASKRSYQFWDTPVPKLGCVNTHGPVEP DKDNIRQEPYTLPGFTWDALDLGDRGVLKELYLLNENYVEDDDNMFRFDYSPEFLLWALRPPGWLPQW HCGVRVSSRKLVGFIAPANIHIYDTEKKMVEINFLCVHKKLRSKRVPVLIREITRRVHLEGIFQAV YTAGVVLKPVGTCRYWHRSLNPRKLIKFSHLSRNMTRMQLYRPLPETPKTAGLRPMETKDIPVVH QLLTRYLKQFHLTPVMSQEEVEHWFYPQENIIDTFVENANGEVTDFLSFYTLPTSTIMNHPTHKSLKAAY SFYNVHTQTPLLDLMSDALVLAKMKGDFVFNALDLMENKTFLEKLFKFGIDGNLQYYLYNWKCPMGAEK VGLVLQ</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	56.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online >](#)

RefSeq: [NP_066565](#)

Locus ID: 4836

UniProt ID: [P30419](#)

RefSeq Size: 4903

Cytogenetics: 17q21.31


RefSeq ORF: 1488

Synonyms: NMT

Summary: Myristate, a rare 14-carbon saturated fatty acid, is cotranslationally attached by an amide linkage to the N-terminal glycine residue of cellular and viral proteins with diverse functions. N-myristoyltransferase (NMT; EC 2.3.1.97) catalyzes the transfer of myristate from CoA to proteins. N-myristoylation appears to be irreversible and is required for full expression of the biologic activities of several N-myristoylated proteins, including the alpha subunit of the signal-transducing guanine nucleotide-binding protein (G protein) GO (GNAO1; MIM 139311) (Duronio et al., 1992 [PubMed 1570339]).[supplied by OMIM, Nov 2008]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified NMT1 protein (Cat# [TP300594]). The protein was produced from HEK293T cells transfected with NMT1 cDNA clone (Cat# [RC200594]) using MegaTran 2.0 (Cat# [TT210002]).