

Product datasheet for **TP300606L**

SAPK4 (MAPK13) (NM_002754) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human mitogen-activated protein kinase 13 (MAPK13), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200606 protein sequence Red =Cloning site Green =Tags(s)

MSLIRKKGFYKQDVNKTAWELPKTYVSPTHVGSAYGVSVCASDAIDKRSGEKVAIKLSRPFQSEIFAKRAY
RELLLLKHMQHENVIGLLDVFTPASSLRNFYDFYLVMPFMQTDLQKIMGMEFSEEKIQYLVIYQMLKGLKY
IHSAGVVHRDLKPGNLAVNEDCELKILDFGLARHADAEMTGYVVTRWYRAPEVILSWMHYNQTVDIWSVG
CIMAEMLTGKTLFKGKDYLDQLTQILKVTGVPGTGFVQKLNDAKAKSYIQSLPQTPRKDFTQLFPRASPQ
AADLLEKMLELDVDRKRLTAAQALHPFFEPFRDPEEETEAQPFDDSLHEKLTVDDEWKQHIIYKEIVNFS
PIARKDSRRRSGMKL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	41.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_002745</u>
Locus ID:	5603



[View online »](#)

UniProt ID: [O15264](#), [A0A024RD04](#)

RefSeq Size: 6348

Cytogenetics: 6p21.31

RefSeq ORF: 1095


Synonyms: MAPK-13; MAPK 13; p38delta; PRKM13; SAPK4

Summary: This gene encodes a member of the mitogen-activated protein (MAP) kinase family. MAP kinases act as an integration point for multiple biochemical signals, and are involved in a wide variety of cellular processes such as proliferation, differentiation, transcription regulation and development. The encoded protein is a p38 MAP kinase and is activated by proinflammatory cytokines and cellular stress. Substrates of the encoded protein include the transcription factor ATF2 and the microtubule dynamics regulator stathmin. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Jul 2012]

Protein Families: Druggable Genome, Protein Kinase

Protein Pathways: Amyotrophic lateral sclerosis (ALS), Epithelial cell signaling in Helicobacter pylori infection, Fc epsilon RI signaling pathway, GnRH signaling pathway, Leukocyte transendothelial migration, MAPK signaling pathway, Neurotrophin signaling pathway, NOD-like receptor signaling pathway, Progesterone-mediated oocyte maturation, RIG-I-like receptor signaling pathway, T cell receptor signaling pathway, Toll-like receptor signaling pathway, VEGF signaling pathway

Product images:



Coomassie blue staining of purified MAPK13 protein (Cat# [TP300606]). The protein was produced from HEK293T cells transfected with MAPK13 cDNA clone (Cat# [RC200606]) using MegaTran 2.0 (Cat# [TT210002]).