

Product datasheet for TP300607L

MSH5 (NM_002441) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human mutS homolog 5 (E. coli) (MSH5), transcript variant 3, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200607 protein sequence Red=Cloning site Green=Tags(s)

MASLGANPRRTPQGPRPGAASSGFSPAPVPGPREAEEEEVEEEEEELAEIHLCVLWNSGYLGIAYYDTS
STIHFMPDAPDHESLKLLQRLVLEINPQSVVTSAKQDENMTRFLGKLASQEHREPKRPEIIFLPSVDFGL
EISKQRLLSGNYSFIPDAMTATEKILFLSSIIPFDCLLTVRALGGLLKFLGRRRIGVELEDYNVSVPILG
FKKFMLTHLVNIDQDTYSVLQIFKSESHPSVYKVASGLKEGLSLFGILNRCHCKWGKLLRLWFTRP
THD LGELSSRLDVIQFFLLPQNLDMAQMLHRLLGHIKKNVPLILKRMKLSHTKVSQVLYKTVYSALGLRDAC
RSLPQSIQLFRDIAQEFSDDLHHIASLIGKVDFEGSLAENRFTVLPNIDPEIDEKRRMLMGLPSFLTEV
ARKELENLDSRIPSCSVIYIPLIGFLLSIPRLPSMVEASDFEINGLDFMFLSEEKLYHSARTKELDALL
GDLHCEIRDQETLLMYQLQCQVLARAVALTRVLDLASRLDVLALASAARDYGYSPRYSPQVLGVRIQN
GRHPLMELCARTFVPNSTECGGDKGRVKVITGPNSSGKSIYKQVGLITFMALVGSFVPAEEAEIGAVDA
IFTRIHSCESISLGLSTFMIDLNQVAVANNATAQSLVLIDEFGKGTNTVDGLALLAAVLRHWLARGPT
CPHIFVATNFLSLVQLQLLPQGPLVQYLTMETCEDGNDLVFFYQVCEGVAKASHASHATAAQGLPKLVA
RGKEVSDLIRSGKPIKPVKDLLKKNQMENCQTLVDKFMKLDLEDPNLDLNVFMSQEVLPAAATSIL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	92.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_002432
Locus ID:	4439
UniProt ID:	O43196 , A0A024RCM1
RefSeq Size:	2945
Cytogenetics:	6p21.33
RefSeq ORF:	2505
Synonyms:	G7; MUTSH5; NG23; POF13
Summary:	This gene encodes a member of the mutS family of proteins that are involved in DNA mismatch repair and meiotic recombination. This protein is similar to a <i>Saccharomyces cerevisiae</i> protein that participates in segregation fidelity and crossing-over events during meiosis. This protein plays a role in promoting ionizing radiation-induced apoptosis. This protein forms hetero-oligomers with another member of this family, mutS homolog 4. Polymorphisms in this gene have been linked to various human diseases, including IgA deficiency, common variable immunodeficiency, and premature ovarian failure. Alternative splicing results multiple transcript variants. Read-through transcription also exists between this gene and the downstream chromosome 6 open reading frame 26 (C6orf26) gene. [provided by RefSeq, Feb 2011]
Protein Families:	Druggable Genome

Product images:



Coomassie blue staining of purified MSH5 protein (Cat# [TP300607]). The protein was produced from HEK293T cells transfected with MSH5 cDNA clone (Cat# [RC200607]) using MegaTran 2.0 (Cat# [TT210002]).