

Product datasheet for TP300639L

Lipoamide Dehydrogenase (DLD) (NM_000108) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human dihydrolipoamide dehydrogenase (DLD), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200639 protein sequence Red=Cloning site Green=Tags(s)

MQSWSRVYCSLAKRGHFNRIHGLQGLSAVPLRITYADQPIDADVTVIGSGPGGYVAAIKAAQLGFKTVCIEKNETLGGTCLNVGCIPSKALLNNSHYHMAHGKDFASRGIEMSEVRLNLDKMMMEQKSTAVKALTGGIAHLFKQNKVVHVNGYGGKITGKNQVTATKADGGTQVIDTKNLIATGSEVTPFPGITIDEDTIVSSTGALSCLKVPEKMMVIGAGVIGVELGSVWQRLGADVTAVEFLGHVGGVGDMEISKNFQIRLQKQGFKFLNLTQVTGATKKSDBGKIDVSI EAASGGKAEVITCDVLLVCIGRRPFTKNLGLLEELGIELDPRGRIPVNTRFQTKIPNIYAIGDVVAGPMLAHKAEDGEIICVEGMAGGAVHIDYNCVPSVIYTHPEVAWVGKSEEQLKEEGIEYKVGKFPFAANSRAKTNADTDGMVKILGQKSTDRVLGAHILGPGAGEMVNEAALALEYGASCEDIARVCHAHPTLSEAFREANLAASFGKSINF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	50.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)


RefSeq:	NP_000099
Locus ID:	1738
UniProt ID:	P09622 , A0A024R713
RefSeq Size:	3613
Cytogenetics:	7q31.1
RefSeq ORF:	1527
Synonyms:	DLDD; DLDH; E3; GCSL; LAD; OGDC-E3; PHE3

Summary: This gene encodes a member of the class-I pyridine nucleotide-disulfide oxidoreductase family. The encoded protein has been identified as a moonlighting protein based on its ability to perform mechanistically distinct functions. In homodimeric form, the encoded protein functions as a dehydrogenase and is found in several multi-enzyme complexes that regulate energy metabolism. However, as a monomer, this protein can function as a protease. Mutations in this gene have been identified in patients with E3-deficient maple syrup urine disease and lipoamide dehydrogenase deficiency. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]

Protein Families: Druggable Genome

Protein Pathways: Citrate cycle (TCA cycle), Glycine, serine and threonine metabolism, Glycolysis / Gluconeogenesis, Metabolic pathways, Pyruvate metabolism, Valine, leucine and isoleucine degradation

Product images:



Coomassie blue staining of purified DLD protein (Cat# [TP300639]). The protein was produced from HEK293T cells transfected with DLD cDNA clone (Cat# [RC200639]) using MegaTran 2.0 (Cat# [TT210002]).